

The PediSpa™

by VITA®
BATH



Model 3036PAS

Specifications

Important Installation Information

The PediSpa Massage Chair is fully adjustable. When fully reclined, there must be an allowance of 17 inches from any obstruction in back of the PediSpa. The warranty does not apply if there is damage due to improper mounting.

Spa Description:

- Construction:** Shell made from Lucite Acrylic, reinforced fiberglass
- Chair:** Comes only in black
Mechanically operated reclining and sliding seat with hand held control
Both right and left arms of chair pivot up and chair fully swivels 360 degrees
Handheld control pad pinpoints vibrating motors to massage specific areas
- Features:** Height adjustable footrest
Extending wash and rinse water spray
Spout is positioned well above water level
Aromatherapy canister & E-Z Clean™ sanitation port
Chromatherapy light
Chrome accents

- Dimensions:** Spa (seat upright): 55" length (140 cm) x 27" width (69 cm) x 54" height (137 cm)
Spa (seat reclined): 74" length (188 cm) x 27" width (69 cm) x 48" height (122 cm)

- Weight:** Spa with Chair: 156 lbs
Crated weight: 180 lbs

- Maximum Capacity:** Five (5) U.S. Gallons

- Air Pump Motor:** 1/2 HP, 110V, 60Hz air therapy pump

- Vibrating Motors:** 12V DC vibrating motor with adapter

- Power Source:** 110V AC/15 Amp, 60 Hz dedicated circuit, GFCI protected, 4 outlet power box included

- Hydro-Therapy Air Jets:** 12 air injectors

- Discharge Pump:** Optional 1/15 HP, 115V, 60Hz

- Electrical Switches:** One (1) On/Off push air activated safety switches for air therapy pump
One (1) On/Off push air activated safety switches for Chromatherapy
One (1) On/Off press switch for discharge pump (Optional)
Handheld remote keypad for chair massage motors

- Installation:** Installation should be made with drain and water connections provided at rear of PediSpa
Hot & Cold water supply lines (1/2") and additional drain line not included
Unit must be attached to the floor—Please consult your local plumber & local codes

Specified model meets or exceeds the following: ANSI Z124.1.2 ASME A112.19.7m

Patent Pending

© 2008 DM Industries, Ltd. All Rights Reserved. Specifications subject to change without notice.

Vita Bath is a registered trademark of DM Industries, Ltd.

2320 NW 147 Street, Miami, Florida 33054

Phone: 800.848.2772x142 Fax: 305-423-2530 • www.vitabath.com



E-Z CLEAN™ Full Flush Sanitizing System by **VITA**® **BATH**

The EZ Clean Full Flush Sanitizing System™ offers you complete peace of mind when it comes to the care of your PediSpa by Vita Bath™. Cleans air lines in seconds!

Sanitizing Instructions:

Step 1

- a. Empty vessel and remove aroma canister.
- b. Close aroma port (replace cap on aroma port).
- c. Turn on blower for 10 seconds. While blower is on, use hand-held spray to rinse vessel.



Step 2

- a. Open aroma port.
- b. Pour 5 to 7 ounces of an EPA registered disinfectant solution such as Ameri-Kleen® Pedicure Spa One-Step.
- c. Close the aroma port.
- d. Turn blower on for 30 seconds.
- e. While blower is on, use hand-held spray to thoroughly rinse the vessel.
- f. Leave blower on for 30 seconds after rinsing the vessel to dry the air lines.
- g. Wipe dry the vessel and you are ready for your next pedicure session.

